

Specification

Selective solid medium for the isolation and primary identification of *Streptococcus agalactiae* from clinical specimens.

Presentation

20 Prepared Plates
90 mm
with: 21 ± 2 ml

Packaging Details

1 box with 2 packs of 10 plates/pack. Single cellophane.

Shelf Life

2 months

Storage

2-14 °C

Composition

Composition (g/l):

Proteose Peptone N.3.....	25.0
Soluble Starch.....	10.0
MOPS.....	11.0
Magnesium sulphate.....	0.2
Metotrexate.....	0.006
Metronidazol.....	0.01
Colistin sulphate.....	0.005
Disodium phosphate.....	5.75
Monosodium phosphate.....	1.5
Crystal Violet.....	0.0002
Agar.....	10.0
Horse Serum.....	50 ml

Description /Technique

Description:

This medium is a modification of the New Granada Medium, which was developed from the GBS Islam Medium. The modifications consist in the suppression of glucose and piruvate, the reduction of the content of starch and an increased volume of horse serum.

Technique:

The sample must be processed and inoculated according to any laboratory methodology or following established standards. If the agar plate has been refrigerated, allow to warm to room temperature before inoculation. Streak sample onto plate and incubate at 35 ± 2°C for 18-24 hours. Anaerobic incubation enhances the orange pigment production and incubation in a CO2 enriched atmosphere increases the colonial size. The specific identification must be done by biochemical and/or serological methods or by genetic technologies.

Note: Do not use *Streptococcus agalactiae* ATCC 13813 for a positive control, as this strain is non-hemolytic in the absence of CAMP factors.

Precautions

For in vitro diagnostic use. Do not reuse. For professional use only.

Do not use the product if it shows evidence of microbial contamination, discoloration, drying, cracking or other signs of deterioration.

Quality control**Physical/Chemical control**

Color : Pale yellow

pH: 7.3 ± 0.2 at 25°C

Microbiological controlSpiral Spreading: Practical range 100 ± 20 CFU. min. 50 CFU (productivity) / 10⁴-10⁶ CFU (selectivity).

Microbiological control according to ISO 11133:2014/A1:2018.

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

5% CO₂ atmosphere. Incubation at 35-37 °C during 24-48 h.**Microorganism***Streptococcus pyogenes* ATCC® 19615*Escherichia coli* ATCC® 25922, WDCM 00013*Streptococcus agalactiae* ATCC® 12386**Growth**

Good - Colourless

Inhibited

Good - Deep Orange

Sterility control

Incubation 48 h at 30-35 °C and 48 h at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

Bibliography

de la ROSA, M., R. VILLAREAL, D. VEGA, C. MIRANDA & A. MARTINEZBROCAL (1983) Granada Médium for detection and identification of group B streptococci. J. Clin Microbiol 18:4: 779-785

de la ROSA, M., M. PEREZ, C. CARAZO, L. PAREJA, J.I. PEIS & F. HERNANDEZ (1992) New Granada Medium for detection and identification of group B streptococci. J. Clin. Microbiol. 30:4:1019-1021

EI AILA, N. A., I. TENCY, G. CLAEYS, B. SAERENS, P. COOLS, H. VERSTRAELEN, M. TEMMERMAN, R. VERHELST & M.

VANEECHOUTTE (2010) Comparison of different sampling techniques and of different culture methods for detection of group B streptococcus carriage in pregnant women. BMC Infectious Diseases 10:285

ISLAM, A.K.M.S. (1977) Rapid recognition of group B haemolytic streptococci. Lancet i:256-257

MERRIT, K., & N. J. JACOBS. (1976) Improved medium for detecting pigment production by group B streptococci. J. Clin Microbiol 4:4:379-380

NOBLE, M.A., J.M. BENT & B. WEST (1983) Detection and identification of group B streptococci by use of pigment production. J. Clin. Pathol. 36:350-352

VOTAVA, M., M. TEIKALOVA, M. DRABKOVA, V. UNZEITIG & I. BRAVENY (2001) Use of GBS media for rapid detection of group B streptococci in vaginal and rectal swabs