

GAS STATION

Test compressed gas microbiological quality



GAS STATION

Test compressed gas microbiological quality

Compatibility :

● 90 or 55mm dishes  100L / min

 Air, nitrogen, argon* and CO₂*
*Flowmeter to be adapted, on request.



From 1 to 6 bars

Easy to use

Protocol is similar to classic air sampling



Tuto available on Youtube

Flowmeter

Visual check of impaction speed



Compact and light :
ideal for limited spaces

Stainless steel
Easy to clean





**20 YEARS EXPERIENCE
IN MICROBIOLOGY EQUIPMENT DESIGN**

**A FULLY COMMITTED TEAM
OF SERVICE ENGINEERS AND MICROBIOLOGISTS**




**MORE THAN 70 INTERNATIONAL DISTRIBUTORS
TRAINED TO THE LATEST A.B.E. TECHNOLOGIES**

**A CONNECTED EQUIPMENT GENERATION
FOR TRACEABILITY AND PRODUCTIVITY**




FOOD SAMPLE PREPARATION



AUTOMATED MEDIA PREPARATION



SAMPLE COLLECTION



BOOK YOUR DEMO !

<https://www.alliance-bio-expertise.com/en/store-locator/>



FOLLOW THE BLOG

alliance-bio-expertise.com



ZI de Courbouton, Le Tremplin - 35480 GUIPRY - FRANCE

E-mail : contact@abioexpertise.com

Phone : +33(0)2.40.51.79.53

alliance-bio-expertise.com